

# System Overview

---

rainbow CT



# System Overview

## 1. Comparative Matrix

	Modalities				Detector	Install Option			Onsite Upgrade
	CBCT	Panoramic	Cephalometric	Stitching	Number of	Floor Stand	Floor Mount	Wall Mount	Cephalometric
2-in-1	O	O	X	O	1	O	X	X	X
3-in-1	O	O	O	O	2				

## 2. Reading Serial Number

① ② ③ ④ ⑤ ⑥  
L DD DX D NN YY

[Sample] L20DXD0221

- ① Month (A~L) – A: January
- ② Production Date
- ③ Model Code (DX: rainbow CT / BC: bright CT)
- ④ Dentium
- ⑤ Production number
- ⑥ Year of production

# Cont'd..





## 3. General View

#	Part Name	Description
①	Vertical	FDD (Focal spot to Detector Distance) * CT, Pan: 600mm (23.6") * Ceph: 1,745mm (68.7")
②	Rotator	
③	Detector (CT+Pan)	CMOS
④	Tube Tank	* Tube Voltage: 60 ~ 100kV $\pm$ 3% * Tube Current: 4 ~ 12mA $\pm$ 3%
⑤	Chin Rest	Automatic (Open-Close)
⑥	Control Panel	Manual
⑦	Lift Column	Min. 1,570mm (61.8") / Max. 2,270mm (89.3")
⑧	Main Base	Floor-stand Only
⑨	Ceph Arm	-
⑩	Status LED	Standby (Green), Ready (Blue), Orange (Exposure)
⑪	Secondary Collimator	-
⑫	Detector (Ceph)	CdTe CMOS






Cont'd..

4. Dimensions

	Dimension (W x D x H)	Option	Figure	
W/ Ceph	1,925 x 1,460 x 2,270mm (75.8 x 57.5 x 89.4")	3-in-1		
W/O Ceph	1,050 x 1,460 x 2,270mm (41.3 x 57.5 x 89.4")	2-in-1		


Cont'd..

5. Crates


	Dimension (W x D x H)	Contents	Net weight	Gross weight	Figure
Box 1	1,260 x 1,130 x 395mm (49.6 x 44.5 x 15.5")	Main Base	44kg (97.0lbs)	60kg (132.3lbs)	
Box 2	1,760 x 720 x 1,408mm (69.3 x 28.3 x 55.4")	Main Body Workstation Others	165kg (363.8lbs)	230kg (507.1lbs)	
Box 3	1,650 x 910 x 910mm (65.0 x 35.8 x 35.8")	Cephalometric Accessories	35kg (77.2lbs)	60kg (132.3lbs)	
TOTAL			244kg (538.0lbs)	350kg (771.6lbs)	-

# Cont'd..

## 6. Workstation










	Specification
Manufacturer	HP
Model	Z4 G4 Workstation
CPU	Intel Xeon W-2223
RAM	16 GB DDR4-2933 Reg ECC RAM
HDD	1TB 3.5 inch HDD
Graphic card	nvidia GeForce RTX-3050 8GB
Ethernet	Intel I210 T1 LAN Card / Serial Port
CD/DVD Drive	9.5mm DVD-Writer 1st ODD
OS	Windows 10 Professional
Figure	

+

	Specification
Manufacturer	HP
Model	EliteDisplay E22 G5
Display size	21.5-in
Native resolution	FHD 1920 x 1080 (16:9)
Viewing angle	178°
Figure	


# Cont'd..

## 7. Software Package

Category	Figure	Function
Software	 rainbow 3D Viewer 1.1.0_rev1721  rainbowBackupRestore_1.0.35.3  rainbowDBServer_1.0.35.3  rainbowImageViewer_win7_package_1.0.35.3  rainbowCT_Ver 1.0.0	* 3D Viewer: CBCT viewer * BackupRestore: DB Table * DBServer: PostgreSQL 9.5 * ImageViewer: 2D viewer (Main viewer) * rainbowCT: Image acquisition
Firmware	 rainbowCT_M1.0.1.2  rainbowCT_S1.0.0.4	* M: Main board * S: Sub board
Geometry	 PanoramicGeometry_2.0.4.9_20201208	All modalities
Engineering	 rainbowCT_Ver.02.15	Device test and update firmware

Cont'd..

8. System Configuration





# Mechanical Specification

## 1. Image Magnification


Modality	FOD(mm)*	FDD(mm)**	ODD(mm)***	Magnification****
CBCT	393	600	207	1.53:1
Panoramic	443		157	1.35:1
Cephalometric	1,525	1,745	220	1.14:1

\* FOD (Focal spot to Object Distance)

\*\* FDD (Focal spot to Detector Distance)



\*\*\* ODD (Object to Detector Distance) = FDD-FOD

\*\*\*\* Magnification = FDD / FOD





# Cont'd..

## 2. X-ray Generator

Item		Specification	Figure
Input power		110~240V, 50/60Hz	-
Rated output power		2.0-2.2 kVA	
High Voltage Generator	kVp	60 ~ 100 kVp $\pm$ 3% (1kV step)	
	mA	4 ~ 12 mA $\pm$ 3% (0.1mA step)	
Tube	Manufacturer	POSKOM	
	Insert	Toshiba	
	Focal spot	0.5mm	
	Target Angle	12 °	
	Anode heat content	30 KJ	
	Inherent filtration	1.0mm Al	

Cont'd..

3. Detector

	Panoramic	CBCT	Cephalometric
Type	CMOS		CdTe CMOS
Model	C12820DK-40		C10502D-43
Pixel size	120μm x 120μm	240μm x240μm	100μm x 100μm
Active area	—		290mm x 6mm
Figure			

# Cont'd..

## 4. Electrical Spec.

	Description
Power supply voltage	110-230 V
Frequency	50 / 60Hz
Power rating	2.2 kVA Max.

This product is designed and manufactured to meet the following standards:

- \* IEC/EN/UL 60601-1, IEC/EN 60601-1-1, IEC/EN 60601-1-2, IEC/EN 60601-1-3, IEC/EN 60601-2-7, IEC/EN 60601-2-28, IEC/EN 60601-2-32, IEC/EN 60601-2-44
- \* ISO 9001, ISO 13485



CE symbol grants the product compliance to the European Directive for Medical Devices 93/42/EEC as amended by 2007/47/EC as a class IIb device.

# Panoramic

## System Default


Scan Type		System Default
Patient size	Male	79kVp / 10mA
	Female	78kVp / 10mA
	Pediatric	71kVp / 10mA
Scan time	-	19 s
Protocol	-	Standard
Pixel size	-	0.12mm

## Examination Programs and Scan Time

Protocol	Patient size	Scan Time(S)	Exposure
Standard	Male	19	79kVp / 10mA
	Female		78kVp / 10mA
	Pediatric	16.5	71kVp / 10mA
Incisor	Male	8	79kVp / 10mA
	Female		78kVp / 10mA
	Pediatric	7	71kVp / 10mA
Bitewing	Male	15	79kVp / 10mA
	Female		78kVp / 10mA
	Pediatric	13	71kVp / 10mA
Sinus	Male	12	90kVp / 10mA
	Female		83kVp / 10mA
	Pediatric		75kVp / 10mA
TMJ	Male	9	79kVp / 10mA
	Female		78kVp / 10mA
	Pediatric	8.5	71kVp / 10mA






# Panoramic

## 1. Standard Mode



### ■ Features



- Preset: Full, Incisor, Bitewing

	Chinrest	Bite	Patient position
Standard			
Edentulous			NA









# Cont'd..

## 2. Sinus Mode





- Features
- N/A



	Chinrest	Bite	Patient position
Standard			 
Edentulous			NA

# Cont'd..

## 3. TMJ Mode



- Features
- N/A

	Chinrest	Bite	Patient position
Standard		X	



# CBCT

## System Default



Scan Type		System Default	Voxel Size (um)	
			Fast Normal	High Res.
Patient size	Male	94kVp / 8mA	300um	200um
	Female			
	Pediatric	94kVp / 7mA		
Scan time	-	20 s		
Protocol	-	Standard		
Pixel size	-	0.24 mm		

## Examination Programs and Scan Time

Mode	Protocol	Voxel (um)	F.O.V (cm)	Scan time (Sec)
Normal	Local	300	5x5	20
	Fast		16x10	10
	Standard			20
	High Res.	200		
Stitching	Fast	300	16x18	20
	Standard			40
	High Res.			

# CBCT

## 1. Local Mode (FOV: 5x5cm)




### ■ Features

- 1) Selection: Maxillary, Occlusion, Mandible
- 2) Region: Single tooth
- 3) Clinical Implication: Single implant, Endo, Perio

	Chinrest	Bite	Patient position
Standard			
Edentulous			NA

# Cont'd..

## 2. Standard Mode (FOV: 16x10cm)





### ■ Features

- 1) Selection: Maxillary, Mandible
- 2) Region: Dual arches back to the 3<sup>rd</sup> molars / Dual TMJs
- 3) Clinical implementation:
  - Surgical guides
  - Simultaneous TMJ diagnosis

	Chinrest	Bite	Patient position
Standard			
Edentulous			NA






# Cont'd..

## 3. Stitching Mode (FOV: 16x18cm)



### ■ Features

- 1) 8.0cm higher FOV than Standard mode
- 2) Region: Dual arches back to the 3<sup>rd</sup> molars/sinus
- 3) Clinical implication : Maxillofacial diagnosis

	Chinrest	Bite	Patient position
Standard			
Edentulous			NA

# Cephalometric

## System Default

Scan Type		System Default
Patient size	Male	95kVp / 10mA
	Female	92kVp / 10mA
	Pediatric	88kVp / 10mA
Scan time	-	8 s
Protocol	-	Lateral ( <b>Half</b> )
Pixel size	-	0.10mm

## Examination Programs and Scan Time


Protocol	Patient size	Scan Time(S)	Exposure
Lateral (Half)	Male	8	95kVp / 10mA
	Female		92kVp / 10mA
	Pediatric		88kVp / 10mA
Lateral (Full)*	Male	19	95kVp / 10mA
	Female		92kVp / 10mA
	Pediatric		88kVp / 10mA
PA / AP	Male	6	95kVp / 10mA
	Female		92kVp / 10mA
	Pediatric		88kVp / 10mA
SMV	Male		95kVp / 10mA
	Female		92kVp / 10mA
	Pediatric		88kVp / 10mA
Carpus	Male		95kVp / 10mA
	Female		92kVp / 10mA
	Pediatric		88kVp / 10mA

\* Full scan mode (19 sec) can be activated upon a customer's request.



# Cephalometric

## 1. Lateral Mode (ROI: 24x29cm)





### ■ Features

- 1) Full scan mode (19 sec) can be activated upon a customer's request.
- 2) Purpose: Check the positional changes of the jaw bone and teeth before and after orthodontic treatment

# Cont'd..

## 2. PA Mode (ROI: 24x29cm)





### ■ Features

- Purpose: Capable of observing soft tissue contours, frontal sinus and nasal cavity

# Cont'd..

## 3. SMV Mode (ROI: 24x29cm)




### ■ Features

- Purpose: Position of the mandibular condyle and observation of the base of the skull



# Cont'd..

## 4. Carpus Mode (ROI: 24x29cm)



### ■ Features

- Purpose: Study changes in the skull or compare teeth by assessing bone age

**Thank you**